Transport Simulation Test Report

Single-Use Shipper with 20 X 20 L CryoVault® Containers

Index


1. Approval Signatures	.03
2. Test Overview	.04
2.1 Objective	.04
2.2 Test Results	.04
2.3 Conclusions	.04
2.4 Recommendations	.04
3. Test Product	.05
4. Test Procedure	.06
4.1 Pre-Test Inspection	.06
4.2 Preconditioning & Packing	.06
4.3 Rotational Corner Drop Test (Method B)	.06
4.4 Rotational Flat Drop Test (Method C)	.07
4.5 Random Vibration Test	.07
4.6 Rotational Corner Drop Test (Method B)	.08
4.7 Rotational Flat Drop Test (Method C)	.08
4.8 Post-Test Inspection	.08
5. Test Results	.09
5.1 Pre-Test Inspection	.09
5.2 Preconditioning & Packing	.09
5.3 Rotational Corner Drop Test (Method B)	.09
5.4 Rotational Flat Drop Test (Method C)	.09
5.5 Random Vibration Test	.10
5.6 Rotational Corner Drop Test (Method B)	.10
5.7 Rotational Flat Drop Test (Method C)	.10
5.8 Post-Test Inspection	.10
Appendix A: Applied Testing Equipment	.12
Appendix B: Random Vibration Test Graphs	.33
Appendix C: History of Change	.39

Approval Signatures

Written: _____ Date: 12 February 2024

Antonio Francisco da Fonseca
Topa Institute, Mechanical Project Engineer

Verified: _____ Date: 12 February 2024

Wout Adriaanse
Topa Institute, Thermal Project Coordinator

Date: 12 February 2024
Report number: T23-1239
Revision: 0
Research Date: 26 January 2024

Customer: Topa Thermal
Torenlaan 10
2215 RW Voorhout
The Netherlands

Test Facility: Topa Thermal
Torenlaan 10
2215 RW Voorhout
The Netherlands

Abstract: This report details the methods, specifications and results for the transport simulation of Single-Use Shipper with 20 x 20 L CryoVault® containers.

2. Test Overview

2.1 Objective

The objective of this test is to investigate the influence of worldwide pallet transport on the Single-Use Shipper with 20 x 20 L CryoVault® containers.

2.2 Test Results

The Single-Use Shipper with 20 x 20 L CryoVault® containers has been tested with the following results:

Test	Results
Preconditioning & Packing	No remarks
Rotational Corner Drop	No remarks
Rotational Flat Drop	No remarks
Random Vibration	No remarks
Rotational Corner Drop	No remarks
Rotational Flat Drop	No remarks

2.3 Conclusions

No visual damage to the outside of the pallet load or on the CryoVault® containers. The Single-Use Shipper with 20 x 20 L CryoVault® containers is suitable for worldwide pallet transport.

2.4 Recommendations

Corrugated cardboard has a lower tolerance for humidity than plastic. Replacing the current CryoVault® handling box from corrugated cardboard to 3 mm corrugated plastic will improve the performance of the handling box.

3. Test Product

The test product consists of:

Insulated Container

Description: One fully conditioned KCVS-DS0001A

Dimensions: 1,570 x 1,100 x 1,610 mm

Weight: 360.0 kg (including dry ice)

The KCVS-DS0001A will contain the following product load:

Product Load

Description: One pallet load containing 20 CryoVault® containers.

Dimensions: 1,200 x 800 x 1,280 mm

Weight: 151.0 kg


Figure 1: Fully packed GTS-1150L4D-DI container

4. Test Procedure

The testing procedure, according ASTM D4169-22, Distribution Cycle 2, consists of the following parts:

4.1 Pre-Test Inspection

Before the start of the test, the bags, components, and the insulated containers are inspected by an engineer from Topa Institute. Any damage or irregularities are noted. All faces of the test pallet are identified according to **Figure 2**, where Face 5 is one of the smallest vertical faces of the pallet.

4.2 Preconditioning & Packing

The CryoVault® containers are conditioned at -80 °C for at least 72 hours prior to packing. Once the CryoVault® containers are conditioned, a pallet load of 20 CryoVault® containers are built and placed inside the container before the start of the test.

4.3 Rotational Corner Drop Test (Method B)

A rotational corner drop test is conducted with a loaded pallet. Four drops are conducted; one drop on each corner of the pallet. The opposite corner of the corner on which the pallet is dropped is supported by a 100 mm thick piece of timber.

The drop height depends on the weight of the pallet. This test is based on ASTM D4169-22 and D6179-07 (Assurance Level II, Method B).


Figure 2: Conditioning setup

Pallet Weight	Assurance Level I	Assurance Level II	Assurance Level III
0 to 226.8 kg	305 mm	229 mm	152 mm
226.8 kg and up	229 mm	152 mm	76 mm

Table: Drop Heights

Figure 3: Rotational Corner Drop setup

Figure 4: Rotational Flat Drop Test

4.4 Rotational Flat Drop Test (Method C)

A rotational flat drop test is conducted with a loaded pallet. Two drops are conducted; one drop on the handling side of the pallet and one drop on the opposite side of the pallet. If the pallet can be handled in all directions, a flat drop on all four sides is conducted. The drop height depends on the weight of the pallet. This test is based on ASTM D4169-22 and D6179-07 (Assurance Level II, Method C).

Pallet Weight	Assurance Level I	Assurance Level II	Assurance Level III
0 to 226.8 kg	305 mm	229 mm	152 mm
226.8 kg and up	229 mm	152 mm	76 mm

Table: Drop Heights

4.5 Random Vibration Test

4.5.1 Truck Transport

In order to simulate vibrations during truck transport, a random vibration test is conducted. Three test profiles from ASTM D4169-22 and D4728-17 are used with a frequency range from 1 Hz to 200 Hz and average loads of 0.40 Grms for 40 minutes, 0.54 Grms for 15 minutes, and 0.70 Grms for 5 minutes. The entire test lasts 60 minutes.

	Truck Low Level	Truck Medium Level	Truck High Level
Average Load	0.40 Grms	0.54 Grms	0.70 Grms
Test Duration	40 minutes	15 minutes	5 minutes

Table: Random Truck Vibration Test Distribution

Figure 5: Random Vibration Test setup

4.5.2 Air Transport

In order to simulate vibrations during air transport, a random vibration test is conducted. Three test profiles from ASTM D4169-22 and D4728-17 are used with a frequency range from 1 Hz to 200 Hz and average loads of 0.16 Grms for 27 minutes, 0.22 Grms for 10 minutes, and 0.29 Grms for 3 minutes. This sequence is repeated three times. The entire test lasts 120 minutes.

	Air Low Level	Air Medium Level	Air High Level
Average load	0.16 Grms	0.22 Grms	0.29 Grms
Test duration	27 minutes (3x)	10 minutes (3x)	3 minutes (3x)

Table: Random Air Vibration Test Distribution

4.6 Rotational Corner Drop Test (Method B)

A rotational corner drop test is conducted with a loaded pallet. Four drops are conducted; one drop on each corner of the pallet. The opposite corner of the corner on which the pallet is dropped is supported by a 100 mm thick piece of timber. The drop height depends on the weight of the pallet. This test is based on ASTM D4169-22 and D6179-07 (Assurance Level II, Method B).

Pallet Weight	Assurance Level I	Assurance Level II	Assurance Level III
0 to 226.8 kg	305 mm	229 mm	152 mm
226.8 kg and up	229 mm	152 mm	76 mm

Table: Drop heights

4.7 Rotational Flat Drop Test (Method C)

A rotational flat drop test is conducted with a loaded pallet. Two drops are conducted; one drop on the handling side of the pallet and one drop on the opposite side of the pallet. If the pallet can be handled in all directions, a flat drop on all four sides is conducted. The drop height depends on the weight of the pallet. This test is based on ASTM D4169-22 and D6179-07 (Assurance Level II, Method C).

Pallet Weight	Assurance Level I	Assurance Level II	Assurance Level III
0 to 226.8 kg	305 mm	229 mm	152 mm
226.8 kg and up	229 mm	152 mm	76 mm

Table: Drop Heights

4.8 Post-Test Inspection

After all tests are performed, the bags, components, and the insulated containers are inspected by an engineer from Topa Institute. Any damage or irregularities are noted.

5. Test results

5.1 Pre-Test Inspection

Before the tests are performed, the outside of the pallet load is visually inspected by an engineer from Topa Institute.

Observations:

No remarks.

5.2 Preconditioning & Packing

The CryoVault® containers are conditioned at -80 °C for at least 72 hours prior to packing.

Observations:

No remarks.

Figures 6 to 9: Packing

5.3 Rotational Corner Drop Test (Method B)

Weight: Product/Packaging Combination	511.0 kg
Drop Height	152 mm
Faces Tested	5, 2, 6, and 4

Table: Rotational Corner Drop Test

Observations:

No remarks.

5.4 Rotational Flat Drop Test (Method C)

Weight: Product/Packaging Combination	511.0 kg
Drop Height	152 mm
Faces Tested	5, 2, 6, and 4

Table: Rotational Corner Drop Test

Observations:

No remarks.

5.5 Random Vibration Test

	Truck Low Level	Truck Medium Level	Truck High Level
Average Load	0.40 Grms	0.54 Grms	0.70 Grms
Test Duration	40 minutes	15 minutes	5 minutes

Table: Random Truck Vibration Test Distribution Results

	Air Low Level	Air Medium Level	Air High Level
Average Load	0.16 Grms	0.22 Grms	0.29 Grms
Test Duration	27 minutes (3x)	10 minutes (3x)	3 minutes (3x)

Table: Random Air Vibration Test Distribution Results

Results graphs for the random vibration tests can be found in Appendix B.

Observations:

No remarks.

5.6 Rotational Corner Drop Test (Method B)

Weight: Product/Packaging Combination	511.0 kg
Drop Height	152 mm
Faces Tested	5, 2, 6 and 4

Table: Rotational Corner Drop Test Results

Observations:

No remarks.

5.7 Rotational Flat Drop Test (Method C)

Weight: Product/Packaging Combination	511.0 kg
Drop Height	152 mm
Faces Tested	5, 2, 6 and 4

Table: Rotational Flat Drop Test Results

Observations:

No remarks.

5.8 Post-Test Inspection

After all tests are performed, the test product is thawed at lab conditions and is visually inspected with the following results:

Observations:

No visual damage to the outside of the pallet load.

No damage is found on the CryoVault® containers.

Figures 10 and 11: Pallet load after testing

The Topa Institute for packaging and distribution is emphasizing its activities on research and the optimization of products and packaging in their performance during their actual lifetime and distribution. Concerning advice and research, the Topa Institute for packaging and distribution is independent from any company, institute, or organization.

Responsibility

The Topa Institute for packaging and distribution or people acting in research activities are not responsible for damage sustained after use of results of research activities.

Appendix A: Applied Testing Equipment

1. Vibration Test System

• Manufacturer:	Lansmont Model 10000 TTV II
• Maximum size:	152 x 152 cm
• Maximum specimen weight:	1,100 kg
• Maximum amplitude (peak to peak):	6.4 cm
• Frequency range:	1 - 300 Hz
• Frequency range at maximum load:	1 - 200 Hz
• Acceleration range:	0 - 8 g
• Maximum acceleration at maximum load:	0.77 g
• Automatic displacement or acceleration control	
• Automatic sweep generator and random vibration facilities (Lansmont TouchTest Vibration system)	
• Accelerometer	
• Calibration due date:	June 2024

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Ref transducer	Manufacturer :	Kistler
	Model	8704B100
	Serial number	2097568
	Customer ID nr	-

Calibration method:

The calibration was performed on site. A reference accelerometer is attached near to the center of the table and the values are compared to the display read out values of the system. Overall values are measured and FFT spectra taken

Uncertainties:

The uncertainties were calculated in accordance with EA-4/02: Expression of Uncertainty of Measurement in Calibration. A coverage factor of 2 sigma ($k=2$) has been applied to the standard uncertainty to express the expanded uncertainty at approximately 95%. Uncertainty 5 %

Traceability:

The measurements have been executed using standards traceable to (inter)national standards. Supporting documentation relative to traceability is on file and is available on request

Environmental conditions:

Air pressure	-	hPa
Temperature	-	°C
Relative humidity	-	%

Date of Receipt	At location in Voorhout
Date of Calibration :	15 juni 2023
Date of Certificate :	15 juni 2023

Ing. P. Smit

Authorized Signatory : Calibration Manager

Visual Inspection

Results

Connectors Ok
Housing mechanics Ok

Reference transducer fixing Wax and cable protectors

The input is generated by applying several known profiles

The values are then compared with a calibrated accelerometer output. Profiles are measured from 1 Hz to 1 kHz as Overall value.

If the significant frequency is known then measured with a FFT analyser from 0 - 500 Hz with a 0,625 Hz bandwidth.

Calibration profile : ASTM D 4169-09 rail level III

Ref value	Displayed value LANSMONT	Deviation	Unit
0,209	0,204	-0,005	g RMS

Calibration profile : ASTM D 4169-09 truck level II

Ref value	Displayed value LANSMONT	Deviation	Unit
0,528	0,521	-0,007	g RMS

Calibration profile : ASTM D 4169-09 air level I

Ref value	Displayed value LANSMONT	Deviation	Unit
1,538	1,510	-0,028	g RMS

Calibration profile : ASTM D 4169-09 air low 2022

Ref value	Displayed value LANSMONT	Deviation	Unit
0,163	0,157	-0,006	g RMS

Calibration profile : ISTA 3 A Pickup and delivery

Ref value	Displayed value LANSMONT	Deviation	Unit
0,460	0,458	-0,002	g RMS

Calibration profile : 1 " loose load test

4 Hz





Ref value 4Hz	Displayed value LANSMONT	Deviation	Unit
0,730	0,820	0,090	g pk

Calibration profile : Sine Sweep 3 - 300 Hz

Frequency in Hz	Ref value in g pk	Displayed value LANSMONT in g pk	Deviation in g pk
5	0,533	0,500	-0,033
10	0,510	0,500	-0,010
20,1	0,501	0,500	-0,001
40,2	0,507	0,500	-0,007
80	0,504	0,500	-0,004
150,1	0,496	0,500	0,004
199,7	0,514	0,500	-0,014
250	0,513	0,500	-0,013
300	0,504	0,500	-0,004

With $g = 9,80665 \text{ m/s}^2$







15.01.15.15.000000000000

Certificate nr.: TC-6013-6

Page 3 of 4

Reproduction of the complete certificate is allowed. Parts of the certificate may only be reproduced with written approval of Sonor Kalibratie Laboratory.

2. Pulley Block

• Manufacturer:	Demag
• Lifting capacity:	1,250 kg
• Hoisting velocity:	4 m/min or 2 m/min
• Year of construction:	2005
• Inspection due date:	April 2024

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3. Pallet Truck With Weighing Scale

• Manufacturer pallet truck:	Toyota BT Lifter LHM200SC
• Manufacturer weighing scale:	Ravas
• Range weighing scale:	0.0 – 2,000.0 Kg ($\pm 0.1\%$)
• Record:	Digital
• Dimensions weighing platforms:	1,190 x 540 mm
• Year of construction:	2019
• Calibration due date:	April 2024

Mettler Toledo B.V.
Franklinstraat 5
4004 JK
TIEL (NL)
0344-638363

Kalibratiecertificaat

Accuracy Calibration Certificate

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4. Ultra Low Freezer #15

• Manufacturer:	NuAire
• Type:	6512E
• Dimensions chamber (width x depth x height):	63 x 60 x 128 cm
• Door dimensions:	63 x 128 cm
• Temperature range:	-85 °C to -40 °C (± 2 °C)
• Programmable microprocessor for temperature:	N/A: constant temperature only
• Calibration due date:	February 2024

Certificate of Calibration

Fluke Nederland B.V.

Certificate Number:	SA01141933	Date of Calibration:	22 Feb 2023
Receive Condition:	IN TOLERANCE	Date of Recalibration:	22 Feb 2024
Return Condition:	IN TOLERANCE	Place of Calibration:	Voorhout
Manufacturer:	NUAIRE	Temperature within:	(23.0 ± 5) °C
Model:	NU-6512E	Humidity within:	(50 ± 30) %rh
Serial Number:	90403726		
Description:	ULTRA LOW FREEZER -40 TO -85 DEGREES		
Procedure:	Manual Procedure		

Customer:	TOPA INSTITUTE VOORHOUT
Customer Asset ID:	FREEZER 15
RMA Number:	606300972

All measurements are traceable to national and/or international standards or have been derived by approved ratio techniques. When possible standards used for this calibration are ISO/IEC 17025 accredited calibrated.

This calibration is performed by a DEKRA certified lab for ISO 9001. This certificate may not be reproduced other than in full. Calibration certificates without signatures, either electronic or handwritten, are not valid.

Issue Date: 22 Feb 2023

Electronically signed

Authorized By
D.W.G. de Visser

Certificate of Calibration

Certificate Number: SA01141933

Remarks

- The calibration status found in this certificate on the top of each results page must be interpreted as:
 - As Found : Data collected before the unit was adjusted and / or repaired
 - As Left : Data collected after the unit has been adjusted and / or repaired
 - Found / Left : Data collected without any adjustment and / or repair performed
- The calibration interval (due date) is the responsibility of the end user.
- According to the European norm 'Operation of electrical installations' NEN-EN 50110-1 release 2013 and the Dutch norm NEN 3140 release 2015 paragraph 5.102.12 through 5.102.16, a safety test is not required. Therefore not performed.
- Temperature conversions (if applicable) are performed according to ISO/IEC 60584:2013 for thermocouples, and ISO/IEC 60751:2022 for resistance temperature devices.

Standards and test-equipment used

Inventory No	Model	Serial No
TS0094	2626-S	C19075
TS0249	TEMPERATURE-PROBE	FTE PT-100 SET

Calibration Data Report

Certificate number: SA01141933

Calibration Date: 22 Feb 23

Supplied value		Location	Lower limit	Measured value	Upper limit	% of Tol.
Temperature (ITS-90) calibration						
setpoint	indicator					
-85 °C	-84 °C	1	-90.0	-84.0 °C	-80.0	1
-85 °C	-84 °C	2	-90.0	-83.8 °C	-80.0	3
Temperature variation over time (30 minutes)						
			-1.0	0.8 °C	1.0	80
setpoint	indicator					
-40 °C	-41 °C	1	-45.0	-40.9 °C	-35.0	3
-40 °C	-41 °C	2	-45.0	-39.6 °C	-35.0	36
Temperature variation over time (30 minutes)						
			-1.0	0.2 °C	1.0	22
Location:						
1			Lower			
2			Upper			

Remarks :

This measurement was done without test object(s) in the chamber.

Customer specification $\pm 2^{\circ}\text{C}$ and $\pm 1^{\circ}\text{C}$ for stability

5. Ultra Low Freezer #20

• Manufacturer:	NuAire
• Type:	6625V30
• Dimensions chamber (width x depth x height):	86 x 60 x 128 cm
• Door dimension:	86 x 128 cm
• Temperature range:	-85 °C to -40 °C (± 2 °C)
• Programmable microprocessor for temperature:	N/A: constant temperature only
• Calibration due date:	February 2024

Certificate of Calibration

Fluke Nederland B.V.

Certificate Number:	SA01140641	Date of Calibration:	20 Feb 2023
Receive Condition:	IN TOLERANCE	Date of Recalibration:	20 Feb 2024
Return Condition:	IN TOLERANCE	Place of Calibration:	Voorhout
Manufacturer:	NEW BRUNSWICK	Temperature within:	(23.0 ± 5) °C
Model:	C660	Humidity within:	(50 ± 30) %rh
Serial Number:	1004-5488-0404		
Description:	ULTRA DEEP FREEZER		
Procedure:	Manual Procedure		

Customer:	TOPA INSTITUTE VOORHOUT
Customer Asset ID:	FREEZER 20
RMA Number:	606300972

All measurements are traceable to national and/or international standards or have been derived by approved ratio techniques. When possible standards used for this calibration are ISO/IEC 17025 accredited calibrated.

This calibration is performed by a DEKRA certified lab for ISO 9001. This certificate may not be reproduced other than in full. Calibration certificates without signatures, either electronic or handwritten, are not valid.

Issue Date: 20 Feb 2023

Electronically signed

Authorized By
D.W.G. de Visser

Fluke Nederland B.V.

E-mail

Telephone

Rev 230217

Brainport Industries Campus 1, 5657 BX, Eindhoven

service.nl@fluke.com

+31 40 267 5300

Page 1 of 3

Certificate of Calibration

Certificate Number: SA01140641

Remarks

- The calibration status found in this certificate on the top of each results page must be interpreted as:
 - As Found : Data collected before the unit was adjusted and / or repaired
 - As Left : Data collected after the unit has been adjusted and / or repaired
 - Found / Left : Data collected without any adjustment and / or repair performed
- The calibration interval (due date) is the responsibility of the end user.
- According to the European norm 'Operation of electrical installations' NEN-EN 50110-1 release 2013 and the Dutch norm NEN 3140 release 2015 paragraph 5.102.12 through 5.102.16, a safety test is not required. Therefore not performed.
- Temperature conversions (if applicable) are performed according to ISO/IEC 60584:2013 for thermocouples, and ISO/IEC 60751:2022 for resistance temperature devices.

Standards and test-equipment used

Inventory No	Model	Serial No
TS0094	2626-S	C19075
TS0249	TEMPERATURE-PROBE	FTE PT-100 SET

Calibration Data Report

Certificate number: SA01140641

Calibration Date: 20 Feb 23

Supplied value	Location	Lower limit	Measured value	Upper limit	% of Tol.	
Temperature (ITS-90) calibration						
setpoint	indicator					
-50.0 °C	-50.0 °C	1	-55.0	-51.0 °C	-45.0	21
-50.0 °C	-50.0 °C	2	-55.0	-51.1 °C	-45.0	22
-50.0 °C	-50.0 °C	3	-55.0	-50.2 °C	-45.0	4
-50.0 °C	-50.0 °C	4	-55.0	-49.7 °C	-45.0	6
-50.0 °C	-50.0 °C	5	-55.0	-50.4 °C	-45.0	8
Temperature variation over time (30 minutes)						
			-1.0	1.0 °C	1.0	99
Temperature (ITS-90) calibration						
setpoint	indicator					
-86.0 °C	-84.0 °C	1	-91.0	-86.2 °C	-81.0	4
-86.0 °C	-84.0 °C	2	-91.0	-86.2 °C	-81.0	4
-86.0 °C	-84.0 °C	3	-91.0	-83.9 °C	-81.0	42
-86.0 °C	-84.0 °C	4	-91.0	-83.5 °C	-81.0	50
-86.0 °C	-84.0 °C	5	-91.0	-84.2 °C	-81.0	36
Temperature variation over time (30 minutes)						
			-1.0	0.1 °C	1.0	10


Location:

- 1 Lower left front
- 2 Upper left rear
- 3 Upper right front
- 4 Lower right rear
- 5 Center


Remarks :

This measurement was done with little test object(s) in the chamber.


Appendix B: Random Vibration Test Graphs


Graph of the Random Vibration Test, Truck Transport Low Level


Graph of the Random Vibration Test, Truck Transport Medium Level


Graph of the Random Vibration Test, Truck Transport High Level


Graph of the Random Vibration Test, Air Transport Low Level, Cycle 1


Graph of the Random Vibration Test, Air Transport Medium Level, Cycle 1


Graph of the Random Vibration Test, Air Transport High Level, Cycle 1


Graph of the Random Vibration Test, Air Transport Low Level, Cycle 2


Graph of the Random Vibration Test, Air Transport Medium Level, Cycle 2


Graph of the Random Vibration Test, Air Transport High Level, Cycle 2

Graph of the Random Vibration Test, Air Transport Low Level, Cycle 3

Graph of the Random Vibration Test, Air Transport Medium Level, Cycle 3

Graph of the Random Vibration Test, Air Transport High Level, Cycle 3

Appendix C: History of Change

0	AF	12 February 2024	All	Create new document