

# UltraCap® H.D.

High Capacity Capsule Filters



# UltraCap® H.D. Capsule Filters

Meissner's UltraCap® H.D. (Heavy Duty) capsule filters are ready-to-use assemblies that offer high flow and throughput with the convenience and cleanliness of a single-use filter assembly.

Designed for processing of medium to large liquid batches, UltraCap® H.D. high capacity capsule filters are optimized for continuous and batch processing in biomanufacturing and for final and prefiltration in pharmaceutical, food and beverage, and microelectronics applications. UltraCap® H.D. filters withstand higher operating pressure and are more robust than conventional UltraCap® capsule filters.

Meissner UltraCap® H.D. capsule filters are optimized for integration into single-use systems such as Meissner's One-Touch® portfolio.

UltraCap® H.D. assemblies are available with a range of Meissner filter media for liquid, gas, and venting applications. They can be specified with a variety of inlet and outlet connections. An optional gauge port facilitates pressure measurement, while an optional filter stand facilitates fast, easy installation.



## Features and Benefits

- Ruggedized polypropylene assembly withstands higher pressures than conventional high capacity capsule filters and resists damage, ensuring reliability and integrity under demanding conditions
- Encapsulated, integral assembly reduces operator contact with filtered liquids
- Extremely low hold-up volume design conserves valuable filtered liquids
- Single-use filter assembly saves installation, setup, cleaning and cleaning validation costs
- Seamlessly integrates into One-Touch® single-use systems or other single-use portfolios
- Removal ratings from 0.04 µm to 99 µm
- Final filtration through prefiltration media options include PVDF, PES, PP and PTFE membranes, as well as PP microfiber, borosilicate glass microfiber, and PP microfiber depth media
- Valved vent port for security and reliability in venting, draining and sampling
- Recessed filter vent/drain on T-style configuration prevents breakage in use
- UltraCap® H.D. filters can be easily configured in series or parallel to maximize design space. 10", 20", 30", 40" and 50" lengths permit fast, easy scale-up
- Can be used with UltraSnap® connectors to configure multiple pre and final capsule filters into a presterilized, ready-to-use assembly
- Available gamma-irradiated for aseptic applications

# Product Specifications

## Materials of Construction

UltraCap® H.D. Housing: Animal component free (ACF), gamma stable polypropylene (PP)

### Filtration Media:

#### Hydrophilic Membranes

SteriLUX®	PVDF (polyvinylidene fluoride)
EverLUX®	PES (polyethersulfone)
STyLUX®	PES (polyethersulfone)

#### Hydrophobic Membranes

Steridyne®	PVDF (polyvinylidene fluoride)
Chemdyne®	PP (polypropylene)
Ultradyne®	PTFE (polytetrafluoroethylene)

#### Microfiber

ALpHA®	PP (polypropylene)
ALpHA® G	PBT (polyester)
Vanguard®	PP (polypropylene)
Protec® RF	GF (borosilicate glass)
Protec® RM	GF (borosilicate glass) + SteriLUX® PVDF membrane
DeltaMax®	PP (polypropylene)
DeltaDepth®	PP (polypropylene)

Support Components: Polypropylene (PP)

Sealing Method: Thermal Bonding

## Cartridge Length (Nominal)

10", 20", 30", 40", or 50"

(25 cm, 50 cm, 75 cm, 100 cm, or 125 cm)

## Max. Pressure & Temperature for Liquids

90 psig @ 32 °F to 100 °F (6.2 bar @ 0 °C to 38 °C)

55 psig @ 140 °F (3.8 bar @ 60 °C)

## Max. Pressure & Temperature for Gases

60 psig @ 32 °F to 100 °F (4.1 bar @ 0 °C to 38 °C)

35 psig @ 140 °F (2.4 bar @ 60 °C)

## Connections

Inlet/Outlet: Sanitary flange, hose barb or Flaretek®

Vent Port: Sanitary valve with hose barb; Sterile Process Design (SPD) vent (inline option only)

Drain Port: Sanitary valve with hose barb; sanitary plug (T-style option only); Sterile Process Design (SPD) vent (inline option only)

Gauge Port: ¾" sanitary flange (T-style option only)

## Sterilization

The UltraCap® H.D. assembly should be autoclaved at a minimum of 121°C for 60 minutes with the vents open to facilitate air removal and the outlet down. UltraCap® H.D. assemblies can be repeatedly autoclaved without loss of integrity. For critical applications, the autoclave cycle should be validated.

UltraCap® H.D. assemblies must not be *in situ* steam sterilized (SIP).

Gamma irradiated models are available.



## Mounting

The UltraCap® H.D. assembly can be mounted and supported on suitably braced, rigid, inline pipe connections. A wall mounting bracket and accessory stand are also available. For applications requiring multiple UltraCap® H.D. capsule filters, Meissner's UltraSnap® filter assembly is recommended. This assembly secures pre and final capsule filters into a single-use filtration system for plug and play use. Contact Meissner for details.





# Configurations



## Ordering Information

UltraCap® H.D. Model	Filter Media - Grade	Retention Rating (µm)	Cartridge Length	Body Style	Inlet/ Outlet	Vent/Drain Ports
<b>CR2</b>	<b>MF</b>	<b>1.2</b>	<b>— 2</b>	<b>T</b>	<b>00</b>	<b>2</b>
<b>CR2</b> = Standard (non-sterile)	<b>Membrane Media</b>	<b>Grade</b>	<b>Retention Rating (µm)</b>	<b>1</b> = 10"	<b>T</b> = T-style	<b>00</b> = 1" sanitary flange
	SteriLUX® PVDF	<b>VTH</b> <sup>1</sup>	0.1, 0.2, 0.4, 0.6	<b>2</b> = 20"	<b>N</b> = Inline	<b>77</b> = 3/4" sanitary flange
<b>GR2</b> = Gamma irradiated		<b>VMH</b> <sup>2</sup>	0.1, 0.2, 0.4, 0.6	<b>3</b> = 30"		<b>02</b> = 1" sanitary flange inlet; 3/8" hose barb outlet
		<b>VLH</b> <sup>3</sup>	0.1, 0.2, 0.4, 0.6	<b>4</b> = 40"		<b>0C</b> = 1" sanitary flange inlet; 1/2" hose barb outlet
	EverLUX® PES	<b>STW</b> <sup>1</sup>	0.2	<b>5</b> = 50"		<b>09</b> = 1" sanitary flange inlet; 9/16" hose barb outlet
		<b>SMH</b> <sup>2</sup>	0.2, 0.4, 0.6			<b>08</b> = 1" sanitary flange inlet; 3/4" hose barb outlet
		<b>SLH</b> <sup>3</sup>	0.2, 0.4, 0.6			<b>0D</b> = 1" sanitary flange inlet; 1" hose barb outlet
		<b>SLW</b> <sup>3</sup>	0.2			<b>22</b> = 3/8" hose barb
	STyLUX® PES	<b>ST</b> <sup>1</sup>	0.04, 0.1, 0.2, 0.4			<b>CC</b> = 1/2" hose barb
		<b>SM</b> <sup>2</sup>	0.04, 0.1, 0.2, 0.4, 0.6			<b>99</b> = 9/16" hose barb
		<b>SL</b> <sup>3</sup>	0.04, 0.1, 0.2, 0.4, 0.6			<b>88</b> = 3/4" hose barb
	Steridyne® PVDF	<b>VTW</b> <sup>1</sup>	0.2			<b>DD</b> = 1" hose barb
		<b>VMV</b> <sup>2</sup>	0.2			<b>AA</b> = 1/2" Flaretek®
						<b>BB</b> = 3/4" Flaretek®
	<i>Gamma-irradiated model not available for media listed below:</i>					
	Chemdyne® PP	<b>PT</b> <sup>1</sup>	0.2			
		<b>PM</b> <sup>2</sup>	0.04, 0.1, 0.2			
	Ultradyn® PTFE	<b>TT</b> <sup>1</sup>	0.2			
		<b>TA</b> <sup>2</sup>	0.2			
		<b>TM</b>	0.05, 0.1, 0.2, 0.4, 1.0, 5.0			
	<b>Microfiber Media</b>	<b>Grade</b>	<b>Retention Rating (µm)</b>			
	ALpHA® PP	<b>MF</b>	0.45, 0.6, 0.8, 1.2, 2.4, 5, 7, 10, 20, 30, 40*, 70*			
	ALpHA® G PBT	<b>MG</b> *	0.6, 1.2, 2.4, 5, 7, 10, 20			
	Vanguard® PP	<b>MN</b>	0.1, 0.2, 0.4, 1, 3, 5, 10, 30, 60, 99 (nominal)			
	Protec® GF	<b>RF</b> *	0.5, 1			
	Protec® GF + PVDF	<b>RM</b> *	0.2, 0.3, 0.5			
	DeltaMax® PP depth	<b>DM</b>	0.5, 1, 3, 5, 10, 20, 40, 70			
	DeltaDepth® PP depth	<b>DD</b>	0.5, 1, 5, 10, 25, 50 (nominal)			
	<i>*ALpHA® G, Protec® RF, Protec® RM, and ALpHA® 40 µm &amp; 70 µm are gamma-irradiatable.</i>					
						<b>T-style</b>
						<b>0</b> = No vent/drain
						<b>1</b> = No vent; 1/4" sanitary drain plug
						<b>2</b> = Sanitary vent; 1/4" sanitary drain plug
						<b>3</b> = Sanitary vent; 3/4" sanitary flange gauge port; 1/4" sanitary drain plug
						<b>4</b> = Sanitary vent; no drain
						<b>5</b> = Sanitary vent; 3/4" sanitary flange gauge port; no drain
						<b>6</b> = No vent/drain; 3/4" sanitary flange gauge port
						<b>A</b> = No vent; sanitary drain valve
						<b>B</b> = Sanitary vent; sanitary drain valve
						<b>C</b> = Sanitary vent; sanitary drain; 3/4" sanitary flange gauge port
						<b>Inline</b>
						<b>0</b> = No vent/drain
						<b>2</b> = Two sanitary vent/drain valves
						<b>4</b> = One sanitary vent or drain valve
						<b>K</b> = Two SPD vent/drain valves at inlet and outlet
						<b>L</b> = One SPD vent at outlet only

### Filter Media Grade Descriptions

<sup>1</sup> T-grade (VTH, STW, ST, VTW, PT, TT)

This absolute, microbially rated filter meets full traceability requirements for the pharmaceutical industry. It is 100% integrity tested during manufacture. Each T-grade filter is shipped with a Certificate of Quality stating exact quality control criteria and test performance results. This is a validatable product to meet the stringent requirements of the pharmaceutical industry.

<sup>2</sup> M-grade (VMH, SMH, SM, VMV, PM, TA)

This absolute, microbially rated filter is 100% integrity tested during manufacture. It is suited for critical applications when regulatory documentation requirements are minimal.

<sup>3</sup> L-grade (VLH, SLH, SLW, SL)

This filter is not 100% integrity tested or flushed during manufacture. It is offered as an economical prefilter or final filter when sterility assurance is not required.

UltraCap® H.D., SteriLUX®, Steridyne®, EverLUX®, STyLUX®, Chemdyne®, Ultradyn®, ALpHA®, Protec®, DeltaMAX®, DeltaDepth® and One-Touch® are registered trademarks of Meissner Filtration Products. Flaretek is a registered trademark of Entegris Inc.

© 2016, 2012, 2007, 2005, 1998 Meissner Filtration Products, Inc. All rights reserved.

CR 3.1



1001 Flynn Road • Camarillo, CA 93012 • USA  
+1.805.388.9911 • +1.805.388.5948 Fax  
www.meissner.com