

## 0.2 µm SM-grade Mini Capsule Filter (CM2/CK2 Model)

### Description

The STyLUX® 0.2 µm is a hydrophilic PES membrane filter compatible with a wide range of liquids. It withstands a wide pH range (1-14) and can be used to remove contaminants from a broad range of pharmaceutical preparations, antibiotics, vaccines, protein solutions, virus suspensions, enzymes, buffers, ophthalmic solutions, reagents, salt solutions, nutrients, serum and blood-based products, biologicals, and wine and beverages. The filter's asymmetric structure provides bacteria and particle removal at high flow rates and low pressure drops. STyLUX® has very low binding characteristics for preservatives commonly used in the pharmaceutical industry and is compatible with most cleaning chemicals, sanitizers and biocides.

### Materials of Construction

All components of the STyLUX® mini capsule are animal component free. These materials are listed for food contact use in the Code of Federal Regulations (CFR), Title 21, as below:

Membrane:	Polyethersulfone	CFR Title 21, 177.2440
Upstream support:	Polypropylene	CFR Title 21, 177.1520
Downstream support:	Polypropylene	CFR Title 21, 177.1520
Outer guard:	Polypropylene	CFR Title 21, 177.1520
Core:	Polypropylene	CFR Title 21, 177.1520
End caps:	Polypropylene	CFR Title 21, 177.1520
Capsule housing:	Polypropylene	CFR Title 21, 177.1520
Sealing method:	Thermal bonding	

**Pore Size** 0.2 µm

**Minimum Bubble Point** 44 psi (3.0 bar), water  
16 psi (1.1 bar), 60% IPA / 40% water

**Bacterial Retention** >10<sup>7</sup> per cm<sup>2</sup> removal of *Brevundimonas diminuta* per ASTM F838

### Operating Characteristics

Operating temperature range: 32 °F to 122 °F (0 °C to 50 °C)  
Maximum operating pressure: 100 psig (6.9 bar)

### Sterilization

Autoclave: 121 °C to 135 °C (15 to 30 psi, 1 to 2 bar), 30 to 60 minutes, ≥ 3 cycles. Water wet membrane prior to autoclaving.  
Gamma irradiation: 25 to 40 kGy once. Do not autoclave irradiated capsules.  
Capsules must not be in-line steam sterilized.


### Biological Safety

STyLUX® filters meet the requirements as specified in the current USP <88> Class VI plastics, physicochemical, oxidizable substances, and USP <87> cytotoxicity tests. Bacterial endotoxin levels in aqueous extracts of STyLUX® filters are less than 0.5 EU/mL, as determined using the *Limulus* amoebocyte lysate (LAL) test. No binders, adhesives or surfactants are used in the construction of STyLUX® filters. Filters comply with European Commission Regulation No 10/2011. STyLUX® CM2/CK2 capsules are TSE/BSE/animal component free (ACF).

### Quality Assurance

STyLUX® filters comply with the Food and Drug Administration Code of Federal Regulations, Title 21, Parts 210 and 211. Product is manufactured and packaged in a cleanroom facility that, through voluntary compliance, meets or exceeds FDA Good Manufacturing Practice Standards. To ensure product reliability, Meissner's Quality Assurance staff continually audits the manufacturing process for conformance to its Quality Management System. Each STyLUX® filter is integrity tested during manufacture and is clearly marked with filter type, lot number and serial number.

## Ordering Guide



Additional information about this filter product is available in the STyLUX® Green Docs document at [www.meissner.com/green-docs](http://www.meissner.com/green-docs).

STyLUX® is a registered trademark of Meissner Filtration Products.  
 © 2020 Meissner Filtration Products, Inc. All rights reserved.