# **AccuFlux**<sup>®</sup>

Integrity Test Instrument



# **MEISSNER**

# **AccuFlux® Integrity Tester**



The AccuFlux<sup>®</sup> integrity tester demonstrates Meissner's commitment to develop smart products that reduce process risk and enhance operator convenience. The AccuFlux<sup>®</sup> is robust and easy to use. It delivers the highest level of accuracy in performing the following filter integrity tests for critical process applications:

- Diffusive flow
   Water intrusion
- Pressure hold

The AccuFlux<sup>®</sup> utilizes three mass flow meters to provide the greatest level of accuracy at all points throughout the standard 0 – 250 mL/ min flow range. These flow meters provide the high level of flow control that allows the diffusive flow test to be performed in less than three minutes on systems containing up to three 40" filter elements without sacrificing accuracy. Bubble point tests use Meissner's TransVector<sup>™</sup> algorithm and a highly accurate pressure sensor to provide test results in 90 seconds.

The AccuFlux<sup>®</sup> features a 15" color multi-touch user interface that is intuitive and easy to use. The AccuFlux<sup>®</sup> system's recipe driven test setup allows users to perform integrity tests seconds after the instrument is turned on. Once powered up, the user logs directly into the test screen, where test recipes are available for selection. From there, a user can review the displayed test data and perform an integrity test with the touch of a button. While the test is being performed, all test settings and data are being stored in up to three locations. Upon completion of a test, the user must acknowledge the result and a test report is generated and saved to a PDF file. Test data can be viewed from the system's History screen. Setup of test recipes is simplified because Meissner filter data is preloaded in the system. The recipes for filters from other manufactures are easily entered manually.

Test data entry and recipe setup is limited to users of a specified permission level. A portion of what makes the AccuFlux<sup>®</sup> system 21 CFR Part 11 compliant is its multi-tiered user administration system. This system requires all users to have a unique user name and password.

The AccuFlux® system's functionality and user experience is enhanced by the following features:

- Preinstalled filter data for Meissner filters
- Recipe driven storage system for storing filter data, system data, and test setups
- Tiered user system, including unique user login accounts
- Multiple selectable save locations, including network locations and USB devices
- Auto-generated, noneditable, PDF test certificate

- Live test data stored to CSV file
- Built-in self checks and tests, including internal/external leak test and source pressure check
- 21 CFR Part 11 Compliant

**Operating System** Windows 10

PC B&R industrial PC 1.91 GHz, quad-core processor 8 GB RAM 128 GB storage

**User Interface** 15.6" color multi-touch display Touch Type: Projected capacitive Resolution: 1920 x 1080 full HD

Connectivity 2.0 USB (3 qty.) 3.0 USB (1 qty.) 24 VDC, 5A power output RJ45, Ethernet WLAN (optional) Bluetooth (optional) Standard Sensor Range Low flow: 0 – 8 mL/min Medium flow: 0 – 40 mL/min High flow: 0 – 250 mL/min Pressure: 0 – 100 PSIG (0-6.9 Bar)

Sensor Accuracy Flow meters:  $\pm 0.5\%$  Rd;  $\pm 0.1\%$  FS Pressure sensor:  $\pm 0.15\%$  FS

Functional Tests Printer test Network test Internal leak test External leak test Source pressure test

Power Requirement 100-240 VAC, 50/60 Hz input power (Internal voltage: 24 VDC) Power: 100 W typical, 200 W peak

Physical Dimensions W x D x H 14.2" x 19.0" x 11.1" (36.0 cm x 48.3 cm x 28.1 cm)

# **Product Specifications**

Pneumatic Requirement

Minimum gas supply pressure: 110 PSIG (7.6 Bar)

Maximum gas supply pressure: 125 PSIG (8.6 Bar)

Gas supply type: Clean dry air

Pneumatic Connections

Air inlet: Color coded, keyed Stäubli coupling

Air outlet: Color coded, keyed Stäubli coupling

Vent inlet: Color coded, keyed Stäubli coupling

Vent outlet: Color coded, keyed Stäubli coupling

Enclosure 304 stainless steel Dust and splash resistant

Environmental Conditions Operating temperature range: 60 °F to 80 °F (16 °C to 27 °C)

Storage temperature range: -4 °F to 122 °F (-20 °C to 50 °C) Maximum humidity: 95% R.H. (No condensation)

The AccuFlux<sup>®</sup> system employs a quad-core, Windows 10 industrial PC with an integrated 15" full color, multi-touch screen. The use of a powerful PC is a must when running today's antivirus software and yields great benefits when integrating the system with a network using standard Windows-based accessories and applications and connecting the system to other network devices and printers.

### The AccuFlux<sup>®</sup> integrity test instrument features:





Integrated calibration functionality

Standard network security and virus protection protocol compatibility

Availability of factory calibration services
 Keyed and color coded process connections
 Unitized manifold design

# AccuFlux<sup>®</sup> Key Features

# Top View Side View

- 15" color multi-touch user interface
- 2 Splash and dust resistant enclosure
- ③ Illuminated power button
- 4 Source air inlet with Stäubli connector
- 5 Air outlet with Stäubli connector
- Exhaust inlet and outlet with Stäubli connector
- High speed USB 3.0

- (8) Universal 100-240V, 50/60 Hz AC power
- 9 Full function Windows 10 PC
- Network & remote accessibility via RJ45
   Ethernet connection
- 1 Auxiliary 5A, 24 VDC power output
- USB accessory port

   (allows connection to barcode scanner, keyboard, mouse, flash drive, wireless network adapter, Bluetooth adapter)

### **Ordering Information**

AccuFlux<sup>®</sup> Model

AX	S234	SE	
AX	<b>S234</b> = Standard flow range 0 - 250 mL/min	SE = Standard splash and dust resistant enclosu	re

AccuFlux<sup>®</sup> is a registered trademark of Meissner Filtration Products, Inc. © 2024, 2023, 2018 Meissner Filtration Products, Inc. All rights reserved. AX 2.4A



(12)

(12)