

Automated Filtration System PUPSIT Skid



Automated Filtration System – PUPSIT Skid

When it comes to sterile filtration of pharma drug substances and drug product, using automation is critical. Worldwide regulatory agencies are promoting automation of repetitive tasks to avoid human errors and guarantee complete batch traceability and accuracy.

Discover your next automated filtration solution with our state-of-the-art Automated Filtration System. Our cutting-edge technology ensures a seamless and rapid Pre-Use Post-Sterilization Filter Integrity Test across a wide range of filter sizes, supporting both single and redundant filter configurations in a single-use closed systems.

Benefits

- cGMP and Annex 1 compliant
- Data collection – 21 CFR Part 21 compliant, batch reports and audit trails linked to customer's SCADA
- Reduced testing time and improved productivity
- Easily configurable to process requirements

Options

1. Available in both single and redundant filtration
2. Level of integration with the SCADA of the production line at the customer's choice
3. Regulation of the filtration flow rate to the needs of your filling equipment
4. Customized Filter and SURF Assemblies and Sampling solutions

Technical Characteristics

Application	• Final Filtration • In-line automated PUPSIT sequence
Filters	• 1 or 2 filter(s) configuration • All standard sizes up to 30" • Customizable slope of filter(s)
Integrity Testing Equipment	Meissner AccuFlux® Technology
Frame	• Stainless Steel 304/316L • Mobile (on wheels)
Flow	Adjustable to the needs of the filling equipment
Dimensions (L x D x H)	Configurable; minimal footprint (850 mm x 1050 mm x 1975 mm)
Electricity	• Power supply: 230V or 120V • IP54
Human Machine Interface	• Machine integration & connection to SCADA/MES • Custom mechanical configuration • Process recipe customization
Data Monitoring	• Flow, pressure, weight • Automatic batch record • Audit/Trail
Interface	Integrated HMI/Touch panel
Certification	• Europe: CE • North-America: UL
Option	Filtration with or without fully automated PUPSIT