QuaDrum®

Lip Clips for stable, secure tubing placement



MEISSNER



The QuaDrum® Lip Clip

The QuaDrum® Lip Clip, which is made of durable stainless-steel, features a clamp with an integrated protective polysulfone insert that gently holds the tubing. It has a clip mount on the clamp's outer edge that attaches to the QuaDrum® rigid outer container to hold the device in place, preventing tubing from slipping during processing.

Designed to enhance functionality and usability, the QuaDrum® Lip Clip helps streamline operations, reduce setup time, and ensure consistent, secure tubing positioning during critical processing steps. The Lip Clip enables quick, tool-free installation while maintaining a firm grip that supports repeatable performance in demanding bioprocessing environments.

Key Features:

- Available in ½" and 1" tubing sizes
- Ergonomic clamping mechanism for quick, tool-free installation
- Integrated mounting feature for ancillary components



Materials of construction

- Tubing Clip Body: 316 Stainless Steel
- Tubing Insert: Polysulfone
- Wing Bolt: 304 Stainless Steel
- Set Screw: 304 Stainless Steel
- Lip Clip: 304 Stainless Steel

Part Numbers:

FASD-001C (½" QuaDrum® Lip Clip)
FASD-002C (1" QuaDrum® Lip Clip)

Quadrum© is a copyright of Meissner Filtration Products, Inc.
© 2026 Meissner Filtration Products, Inc. All rights reserved.
QLC 1.0

Meissner Headquarters
1001 Flynn Road • Camarillo, CA 93012 • USA
+1.805.388.9911 • www.meissner.com

MEISSNER